

## LMD-X2700MD

27-inch 4K 2D LCD medical monitor



### Overview

#### **4K image quality with intuitive operation and flexible installation**

The LMD-X2700MD is a 27-inch 4K LCD medical-grade surgical monitor intended for use in hospital operating rooms, surgical centres, clinics, doctors' offices and similar medical environments. Suitable for use in minimally invasive surgical procedures, it displays bright, high quality 2D colour images with true 4K resolution from endoscopic/laparoscopic cameras and other compatible medical imaging systems.

4K picture quality is enhanced further by HDR gamma support and a wide colour gamut for clearer viewing of fine details and subtle tonal differences. Thanks to Sony's unique upscaling technology, HD/SD resolution signals are automatically upscaled to a 4K view with four times the pixel count of Full HD. HD system customers who invest in 4K displays today, can not only achieve improved image quality but also become ready for eventually upgrading to a 4K system.

Operation is simplified by the monitor's intuitive control panel and LED navigation. There's a wide range of input signal options including HDMI and Display Port, and a choice of picture modes for greater flexibility in today's surgical applications and integrated OR environments.

The compact and easy to clean ergonomic design of the LMD-X2700MD is ideal for integration in modern clinical environments, while installation is simplified with a choice of AC or DC powering options and neat cable management.

## Features

### **Bright, high quality 4K images**

The LMD-X2700MD displays finely-detailed 4K Ultra HD images with a high brightness for superior visibility in the operating room.

### **4K upscaling**

Lower-resolution input signals are upscaled to 4K resolution, offering a sharp, natural view with four times the pixel count of Full HD.

### **HDR Gamma support with HLG**

When connected to a HLG (Hybrid Log-Gamma) input signal from a HLG compatible imaging system, the LMD-X2700MD can display HDR images, accurately reproducing an extended range of brightness levels without clipping white highlights or losing dark shadow details.

### **Wide colour gamut**

The monitor is compliant with BT.2020 colour standards, accurately displaying a significantly wider range of colours than the BT.709 space.

### **Intuitive operation**

The monitor's intuitive control panel with LED navigation allows easy operation, even in dark environments. Three custom buttons can be assigned to frequently used functions.

### **Compact and ergonomic with easy-clean design**

The slim, compact LMD-X2700MD features an ergonomic, easy to hold design that allows simple adjustment of monitor position by hand. Hygiene is simplified by the edge-to-edge flush surface

design that allows liquids and gels to be wiped easily off the LCD panel and controls.

### **Wide range of inputs**

The monitor accepts a wide range of signal input including Display Port, HDMI, DVI without the need for additional input converters.

### **Auto Input Select**

If signals from one input port are accidentally interrupted, the monitor's input switches automatically to the second port as a fail-safe measure.

### **Range of display modes**

A variety of display modes including Mirror Image, Flip image, Picture-in-Picture and Picture-out-Picture are instantly available via push-button selection.

### **Greater installation flexibility**

The monitor can be powered by either direct AC or an external AC adaptor. All signal and power connectors face downwards, simplifying neat installation and cabling.

### **VESA Mounting**

The VESA-mounting standard (100 x 100 mm) simplifies integration in a variety of medical environments.

### **Compliance with medical standards**

This product is distributed to the US and EU as a medical device and satisfies product safety standards (e.g. IEC 60601-1). For more details, please contact your nearest Sony sales office or an authorised dealer.

|                                     |  |
|-------------------------------------|--|
| Panel                               | a-Si TFT Active Matrix LCD   |
| Picture Size (Diagonal)             | 684 mm (26.93 inches)  |
| Effective Picture Size (H x V)      | 596.2 x 335.3, 684.0 mm (w/h, dia)<br>23 1/2 x 13 1/4, 27 inches   |
| Pixel pitch                         | 0.15525 x 0.15525 mm   |
| Resolution                          | H 3,840 dots, V 2,160 lines  |
| Aspect                              | 16:9   |
| Pixel Efficiency                    | 99,99 %  |
| Backlight                           | LED  |
| Luminance (Panel Specification)     | 800 cd/m2 (typical)  |
| Contrast Ratio                      | 1000:1   |
| Colors                              | Approx.1.07 billion colors   |
| Viewing Angle (Panel Specification) | 89°/89°/89°/89° (typical)<br>(up/down/left/right, contrast > 10:1) |
| Gamma                               | 1.8, 2.0, 2.2, 2.4, 2.6, DICOM, HLG                                |

## Input

HDMI connector (x1) HDMI 2.0b,

|                     |  |
|---------------------|--|
| HDMI Input          | HDCP1.4 correspondence   |
| DVI-D Input         | DVI-D connector (x1) TMDS single link                              |
| Display Port        | Display Port connector (x2) 1.2, SST, HDCP1.3 correspondence       |
| Serial Remote (LAN) | D-sub 9-pin (RS-232C) (x1) RJ-45 modular connector (ETHERNET) (x1) |
| Remote              | Stereo mini jack (x1)  |
| AC Input            | AC input connector (x1) 100 V to 240 V, 50/60 Hz                   |
| DC Input            | DC input connector (x1) DC 26 V                                    |

## Output

|                     |                             |
|---------------------|-----------------------------|
| DVI-D Output        | DVI-D connector (x1)        |
| Display Port Output | Display Port connector (x1) |
| DC 5 V Output       | 5 V Output (x1) up to 2A    |
| DC 12 V Output      | 12 V Output (x1) up to 2.5A |

## General

AC IN: 100 V - 240 V, 50/60

|                                      |  |
|--------------------------------------|--|
| Power Requirements                   | Hz, 1.5 A – 0.7 A<br>DC IN: 26 V, 5.2 A (Supplied from AC adaptor)   |
| Power Consumption                    | Approx. 149 W (max.)   |
| Operating Temperature                | 0°C to 35°C<br>32°F to 95°F  |
| Operating Humidity                   | 30% to 85% (no condensation allowed)   |
| Storage/Transport Temperature        | -20°C to +60°C<br>-4°F to +140°F   |
| Storage/Transport Humidity           | 20% to 90%   |
| Operating/Storage/Transport Pressure | 700 hPa to 1060 hPa  |
| Dimensions (W x H x D) *             | 658.8 x 426.8 x 80 mm<br>658.8 x 508.6 x 319 mm (with SU-600 optional stand)<br>26 x 16 7/8 x 3 1/4 inches<br>26 x 20 1/8 x 12 5/8 inches (with SU-600 optional stand) |
| Mass                                 | Approx. 8.8 kg (when the optional stand is not installed)<br>Approx. 19 lb 6.4 oz (when  |

|                      |   |
|----------------------|---|
|                      | the optional stand is not installed)  |
| Mounting             | 100 x 100 mm VESA   |
| Supplied Accessories | AC plug holder (x2)<br>Before Using This Unit (x1)<br>CD-ROM (including the Instructions for Use) (x1)<br>Service Contact List (x1) |
| Optional Accessories | AC adaptor AC-300MD<br>Monitor stand SU-600MD<br>Foot switch FS-24  |
| *1                   | The values for dimensions are approximate.  |

## Related products



### AC-300MD

Provides power for Sony LCD medical monitors



### NUCLEUS

The smart digital imaging platform for medical environments



### LMD-X2710MD

27-inch 4K 2D LCD medical monitor with 12G-SDI

## Gallery

