### LMD-XH550MT

55-inch 4K 3D/2D LCD medical monitor



### Specifications

Picture Performance		
Panel	TFT Active Matrix LCD	
Picture Size (Diagonal)	1387.8 mm (54 3/4 inches)	
Effective Picture Size (H x V)	1209.6 x 680.4mm (47 5/8 x 26 7/8 inches)	
Pixel pitch	0.315 x 0.315 mm (0.0124 x 0.0124 inches)	
Resolution (H x V)	3,840 x 2,160 pixels	
Aspect	16:9	
Pixel Efficiency	99.99%	
Backlight	LED	
Panel Technology	LCD with IPS	
Luminance (Panel Specification)	550cd/m2 (Typical) 1500cd/m2 (Typical, Peak)	

Contrast Ratio	1,000,000:1 (Typical)
Colors	10 bit colors (1,073,741,824)
Panel Frame Rate	50/60 Hz
Viewing Angle (Panel Specification)	89°/89°/89°/89° (Typical) (up/down/left/right, contrast > 10:1)
Vertical Viewing Angle (3D Mode)	37°at a viewing distance more than 1,000mm, crosstalk ratio less than 7% (Typical)
Gamma	1.8, 2.0, 2.2, 2.4, 2.6, DICOM, Highlight, HLG
Input	
HDMI Input	HDMI (x1) (HDCP 2.3 correspondence)
DVI-D Input	DVI-D (x1) (TMDS single link, HDCP 1 4 correspondence)

1.4 correspondence)

SDI Input

Display Port

BNC (x2) 12G/3G/HD-SDI

BNC (x1) 3G/HD/SD-SDI

Display Port (x1) (SST HDCP 1.3 correspondence)	
D-sub 9-pin (RS-232C) (x1), RJ-45	

Serial Remote (LAN)	(x1) (Ethernet, 10BASE-T/100BASE- TX)
Remote	Stereo mini jack (x1)
AC input	AC input connector (x1) AC 100-240 V, 50/60 Hz
DC Input	_
Output	
SDI Output	BNC (x2) 12G/3G/HD-SDI BNC (x1) 3G/HD/SD-SDI
CLONE output	BNC (x1) 12G/3G-SDI
DC 5 V / 12 V Output	5 V Output (x1) 2.0 A max 12 V Output (x1) 2.5 A max
General	
Power Requirements	AC IN: 100 V - 240 V, 50/60 Hz, 3.1 A - 1.1 A
Power Consumption	Approx. 290 W (max.)
Operating Temperature	0°C to 40°C (32°F to 104°F)
Operating Humidity	30% to 85% (no condensation)

Storage/Transport Temperature	-20°C to +60°C (-4°F to +140°F)
Storage/Transport Humidity	20% to 90%
Operating/Storage/ Transport Pressure	700 hPa to 1060 hPa
Dimensions (W x H x D)	1249.6 x 763.3 x 92.8 mm (49 1/4 x 30 1/8 x 3 3/4 inches)
Mass	Approx. 27 kg (Approx. 59 lb 8.4 oz)
Mounting (W x H)	200 x 200 mm 400 x 200 mm
IP Code	IP45 (front)/IP32 (others)

#### Supplied Accessories

AC power cord (1)
Plug holder for the supplied AC
power cord (2)
Before Using This Unit (1)
CD-ROM (including the
Instructions for Use) (1)
Service Contact List (1)
M6 x 12 mm Screw (4)
3D Eye Shield Kit: CFV-E30SK (1)
Instructions for Use of the 3D Eye

Shield Kit (1)

#### **Optional Accessories**

3D Eye Shield Kit: CFV-E30SK 2D Eye Shield Kit: CFV-E20SK 3D Eye Shield: CFV-E30D 2D Eye Shield: CFV-E20D Shield Cover: CFV-50SC Shield Frame: CFV-B100 Foot Switch FS-24

### Related products



LMD-XH320MT

32-inch 4K 3D/2D LCD medical monitor



LMD-XH550MD

55-inch 4K 2D surgical monitor



4000MT

4K 2D/3D medical


recorder

HVO- M

MCC-1000MD

00000

Two-piece Full HD surgical video camera



© 2004 - 2024 Sony Corporation. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features and specifications are subject to change without notice. The values for mass and dimension are approximate. All trademarks are the property of their respective owners.