

CFV-E30SK

3D Eye Shield Kit



Overview

This 3D eye shield starter kit includes one frame (CFV-B100) and 3x disposable 3D eye shields.

The disposable eye shields are intended for use by surgeons and nursing staff in the Operation Room viewing images on the Sony Medical Monitor.

The 3D circular polarising disposable eye shield fits into the shield frame that is light and comfortable to wear for extended periods. Wearing the eye shield and frame, healthcare professionals are protected from blood splashes, body fluids typically encountered in the Operation Room during clinical procedures. The 3D eye shield is compatible with the BKM-30GM and BKM-31GM (can use instead of 30/31GM).

Specifications

Specifications

| | |
|-----------------------|---------------|
| Operating Temperature | 0 °C to 40 °C |
|-----------------------|---------------|

| | |
|--------------------|------------|
| Operating Humidity | 30% to 85% |
|--------------------|------------|

| | |
|-------------------------------|------------------|
| Storage/Transport Temperature | -20 °C to +60 °C |
|-------------------------------|------------------|

| | |
|--------------------------------------|---|
| Storage/Transport Humidity | 20% to 90% |
| Operating/Storage/Transport Pressure | 700 hPa to 1060 hPa |
| Dimensions (W x H x D) *1 | 160 x 45 mm (Shield Frame) (1) |
| | 235 x 77 mm (3D shield) (3) |
| | 6 3/8 × 1 13/16 inches (Shield Frame) (1) |
| | 9 3/8 × 3 1/8 inches (3D shield) (3) |
| Mass | Approx.16 g (Shield Frame) (1) |
| | Approx 7 g (3D shield) (3) |
| | Approx. 0.56 oz (Shield Frame) (1) |
| | Approx 0.25 oz (3D shield) (3) |
| Supplied Accessories | Operating Instructions (1) Attachment Procedure (1) |
| Safety Regulation | This product is distributed to the EU as a medical device and a personal protective equipment. It complies with the related EU regulations. |

Notes

*1

The values for dimensions are approximate.

Related products



CFV-50SC

Shield Cover

Gallery

