

LMD-1530MD

15-inch WXGA 2D LCD medical monitor



Overview

The LMD-1530MD is a widescreen HD medical-certified 15-inch* monitor, with a WXGA panel providing 1280 x 768 pixel resolution. The outstanding image quality and lightweight compact design allows for easy integration into a variety of applications, both inside and outside the operating room. Incorporating an advanced LCD panel with wide viewing angle, it reproduces bright and high-contrast images even in dark areas of the picture. It is equipped with a 10-bit signal processor to produce lifelike images. Smooth gradation and accurate colour representation ensures it is true to the original signals. The monitor accepts a wide variety of SD signal inputs, plus HD input via either HDMI or an analogue component. Its versatile inputs, convenient features, and space-saving design make it an ideal, cost-effective choice for today's demanding medical environments.

* 15.3-inch measured diagonally.

Features

15" Wide Screen LCD Panel

The LMD-1530W is equipped with the latest generation of LCD panel providing high contrast pictures, detailed images with an excellent viewing angle. It offers same picture height as previous Sony CRT 14-inch which makes this product ideal for CRT replacement for entry level applications.

10-bit digital video signal processor for excellent picture reproduction

Incorporating a WXGA-type LCD with 1280 x 768 pixel resolution and a wide viewing angle, the LMD-1530MD offers clear and precise images, essential for today's demanding medical applications. The monitor's 10-bit digital video signal processor reproduces life-like images with accurate color reproduction and a consistent grayscale for static images as well as for moving content.

Anti-reflection (AR) coated protection panel

The LMD-1530MD features a robust AR-coated protection panel to protect the surface of the display and to minimize the damage caused by scratches. It reduces reflections from ambient light even in brightly illuminated operating rooms and clinical areas.

Input versatility

The LMD-1530MD comes equipped with a wide variety of SD inputs including analog composite, Y/C (S-Video), 525i/625i component, and RGB. In addition, it accepts HD signals via HDMI interface or analog component connectors.

Scan mode

The LMD-1530MD has different scan modes for a variety of different image sizes. For ultimate flexibility, users can easily select the most appropriate scan mode for their application, including Full Scan mode in which an image is displayed fully in both vertical and horizontal directions.

Large number of operational convenience features

Five gamma settings

Colour temperature settings

On-screen display (OSD) in seven languages

Built-in tally lamp to identify when images are captured or recorded

Key inhibit function

Parallel remote control interface
 Built-in monitor stand
 VESA 100 mounting (100 x 100 mm spacing)

Compliance with medical standards

This product is distributed to the US and EU as a medical device and satisfies product safety standards (e.g. IEC 60601-1). For more details, please contact your nearest Sony sales office or an authorized dealer.

Specifications

Picture performance

| | |
|----------------|---|
| LCD panel | a-Si TFT Active Matrix |
| Picture size | 334 × 200, 390 mm (W/H, Diagonal) (13 1/4 × 7 7/8, 15 3/8 inches) |
| Resolution | 1280 × 768 dots (WXGA) |
| Viewing angle | 88°/ 88°/ 88°/ 88° (typical) |
| Scan Normal | 0%, Over 5% |
| Aspect | 15:9 |
| Display colour | 16,770,000 |

LINE A input connectors

| | |
|-----------|--------------------|
| Y/C input | Mini DIN 4-pin (1) |
|-----------|--------------------|

| | |
|--|---|
| VIDEO input | BNC (1), 1 Vp-p \pm 3 dB, negative synchronization |
| AUDIO input | Phono jack (1), -5 dBu 47 kilo ohms or higher |
| RGB/COMPONENT input connectors | BNC (3) |
| RGB input | 0.7 Vp-p \pm 3 dB, (Sync On Green, 0.3 Vp-p negative sync.) |
| Component input | 0.7 Vp-p \pm 3 dB, (Sync on Y, 0.3 Vp-p negative sync.) |
| AUDIO input | Phono jack (1), -5 dBu 47 kilo ohms or higher |
| OPTION IN connector | D-sub 9-pin (1), female |
| OPTION AUDIO IN connector | Phono jack (1), -5 dBu 47 kilo ohms or higher |
| External synchronized input connector | BNC (1), 0.3 to 4 Vp-p negative polarity binary |
| HDMI IN connector | HDMI (1) |
| PARALLEL REMOTE input connector | |
| Parallel remote | Modular connector 8-pin (1) |

LINE A output connectors

| | |
|--------------|--|
| Y/C output | Mini DIN 4-pin (1), Loop-through, with 75 ohms automatic terminal function |
| VIDEO output | BNC (1), Loop-through, with 75 ohms automatic terminal function |
| AUDIO output | Phono jack (1), Loop-through |

RGB/COMPONENT output connectors

| | |
|----------------------|---|
| RGB/Component output | BNC (3), Loop-through, with 75 ohms automatic terminal function |
| AUDIO output | Phono jack (1), Loop-through |

External synchronized output

| | |
|-------------------------|---|
| Connector | BNC (1), Loop-through, with 75 ohms automatic terminal function |
| Built-in speaker output | 0.5 W (mono) |

General

| | |
|-------------------|-----------------------------|
| Power | AC 100 V to 240 V, 50/60 Hz |
| Power consumption | |

Maximum: approx. 40 W, 0.6 A to 0.4 A

Dimensions (w/h/d)
Approx. 372 × 336 × 264 mm (14 3/4 × 13 1/4 × 10 1/2 inches) (not including the projection parts)

Mass Approx. 6.2 kg (13 lb 11 oz)

Operating conditions

Temperature 0 °C to 35 °C (32 °F to 95 °F)

Recommended temperature 20 °C to 30 °C (68 °F to 86 °F)

Humidity 30 % to 85 % (no condensation)

Pressure 700 hPa to 1060 hPa

Storage and transport conditions

Temperature -20 °C to +60 °C (-4 °F to +140 °F)

Humidity 0 % to 90 %

Pressure 700 hPa to 1060 hPa

Supplied Accessories

AC power cord 1

| | |
|-----------------------------------|---|
| AC plug holder | 2 |
| Instruction for Use | 1 |
| CD-ROM | 1 |
| Using the CD-ROM Manual | 1 |
| Quick reference | 1 |
| When you First Use the Monitor | 1 |
| Sales Companied Guide | 1 |
| Warranty book | 1 |

Gallery

