

## CFV-E20D

### 2D Eye Shield



#### Overview

This pack contains 15x disposable 2D eye shields. Each shield is intended for use with the CFV-B100 eye shield frame (available separately).

The disposable eye shields are intended for use by clinical staff supporting the surgeon in the Operation Room who may experience discomfort or motion sickness when viewing stereoscopic images on the Sony Medical Monitor. The 2D (left polarising) filter allows Operation Room staff to view monitor images in 2D, offering more wearer comfort without '3D sickness'.

The eye shield fits into the shield frame that is light and comfortable to wear for extended periods. Wearing the eye shield and frame, healthcare professionals are protected from blood splashes, body fluids typically encountered in the Operation Room during clinical procedures.

The 2D eye shield is compatible with the BKM-30GM and BKM-31GM (can use instead of 30/31GM).

#### Specifications

#### Specifications

Operating Temperature	0 °C to 40 °C
Operating Humidity	30% to 85%
Storage/Transport Temperature	-20 °C to +60 °C
Storage/Transport Humidity	20% to 90%
Operating/Storage/Transport Pressure	700 hPa to 1060 hPa
Dimensions (W x H x D) *1	235 x 77 mm (2D shield) (15) 9 3/8 × 3 1/8 inched (2D shield) (15)
Mass	Approx 7 g (2D shield) (15) Approx 0.25 oz (2D shield) (15)
Supplied Accessories	Operating Instructions (1) Attachment Procedure (1)
Safety Regulation	This product is distributed to the EU as a medical device and a personal protective equipment. It complies with the related EU regulations.

## Notes

\*1

The values for dimensions are approximate.

---

## Gallery

