

## CFV-E20SK

#### 2D Eye Shield Kit



#### Overview

This 2D eye shield starter kit includes one frame (CFV-B100) and 3x disposable 2D eye shields.

The disposable eye shields are intended for use by clinical staff supporting the surgeon in the Operation Room who may experience discomfort or motion sickness when viewing stereoscopic images on the Sony Medical Monitor. The 2D (left polarising) filter allows Operation Room staff to view monitor images in 2D, offering more wearer comfort without '3D sickness'.

The eye shield fits into the shield frame that is light and comfortable to wear for extended periods. Wearing the eye shield and frame, healthcare professionals are protected from blood splashes, body fluids and typically encountered in the Operation Room during clinical procedures.

The 2D eye shield is compatible with the BKM-30GM and BKM-31GM (can use instead of 30/31GM).

C	- : C: -	: _	
Spe	CITIC	atio	ns

Operating Temperature	0°C to 40°C

# SONY

Operating Humidity	30% to 85%
Storage/Transport Temperature	-20 °C to +60 °C
Storage/Transport Humidity	20% to 90%
Operating/Storage/Transport Pressure	700 hPa to 1060 hPa
Dimensions (W x H x D) *1	160 x 45 mm (Shield Frame) (1) 235 x 77 mm (2D shield) (3) 6 3/8 × 1 13/16 inches (Shield Frame) (1) 9 3/8 × 3 1/8 inched (2D shield) (3)
Mass	Approx.16 g (Shield Frame) (1) Approx 7 g (2D shield) (3) Approx. 0.56 oz (Shield Frame) (1) Approx 0.25 oz (2D shield) (3)
Supplied Accessories	Operating Instructions (1) Attachment Procedure (1)
	This product is distributed to the EU as a medical

# SONY

Safety Regulation	device and a personal
	protective equipment. It
	complies with the related
	EU regulations.

Notes	
*1	The values for dimensions are approximate.

# SONY

### Gallery



