

CFV-B100

Shield Frame



Overview

This wearable shield frame is intended for use with 3D/2D eye shield and accessories (eye shield is not included with this frame and is available separately).

The disposable eye shields are intended for use by surgeons and nursing staff in the Operation Room viewing images on the Sony Medical Monitor.

The eye shield fits into the CFV-B100 shield frame that is light and comfortable to wear for extended periods. Wearing the eye shield and frame, healthcare professionals are protected from blood splashes, body fluids typically encountered in the Operation Room during clinical procedures.

Specifications

Specifications

| | |
|-----------------------|---------------|
| Operating Temperature | 0 °C to 40 °C |
|-----------------------|---------------|

| | |
|--------------------|------------|
| Operating Humidity | 30% to 85% |
|--------------------|------------|

| | |
|-------------------------------|------------------|
| Storage/Transport Temperature | -20 °C to +60 °C |
|-------------------------------|------------------|

| | |
|--------------------------------------|---|
| Storage/Transport Humidity | 20% to 90% |
| Operating/Storage/Transport Pressure | 700 hPa to 1060 hPa |
| Dimensions (W x H x D) *1 | 160 x 45 mm 6 3/8 × 1 13/16 inches |
| Mass | Approx. 16 g Approx. 0.56 oz |
| Supplied Accessories | Operating Instructions (1) Attachment Procedure (1) |
| Safety Regulation | This product is distributed to the EU as a medical device and a personal protective equipment. It complies with the related EU regulations. |

Notes

*1 The values for dimensions are approximate.

Gallery



