

CCMC-T05

5 metre camera cable for PMW-10MD



Overview

Extends flexibility of medical camera

The CCMC-T05 camera cable allows you to extend the PMW-10MD medical two-piece camera. The PMW-10MD has a compact camera head unit, that can be attached easily either to a surgical microscope, a surgical boom arm, or to a ceiling, and a separated camera control unit. The 5 m cable increases the distance between the camera head and the camera control unit.

MDD approved

The camera cable meets MDD approval in the EU.

Features

Extends distance between camera head and camera control unit

The 5 m camera cable is fully MMD approved in the EU and extends the flexibility of the PMW-10MD in surgical conditions.

Gallery

