

CCMC-SA15

15m standard camera cable for medical video cameras



Overview

The CCMC-SA15 is an optional camera cable that connects the head unit and control unit (CCU) of Sony medical video cameras over distances of up to 15m. Compatible models are the MCC-500MD HD video camera and the MCC- 1000MD Full HD video camera.

Gallery



© 2004 - 2026 Sony Corporation. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features and specifications are subject to change without notice. The values for mass and dimension are approximate. All trademarks are the property of their respective owners.