

Betacam SP

For over twenty years, broadcasters and producers have put their faith in the exceptional quality and performance of Betacam SP tapes.



Overview

Specifications

General Characteristics	
Master Carton Quantity (piece)	BCT-90MLA: 10 BCT-60MLA: 10 BCT-30MLA: 10 BCT-30MA: 50 BCT-20MA: 50 BCT-10MA: 50 BCT-5MA: 50
Master Carton Weight (kg)	BCT-90MLA: 8.4 BCT-60MLA: 7.8 BCT-30MLA: 7.1 BCT-30MA: 16.8 BCT-20MA: 15.8 BCT-10MA: 15.1 BCT-5MA: 14.6
Playing time (min)	BCT-90MLA: 90 BCT-60MLA: 60 BCT-30MLA: 30 BCT-30MA: 30 BCT-20MA: 20 BCT-10MA: 10 BCT-5MA: 5
Tape length (m)	BCT-90MLA: 670 BCT-60MLA: 457 BCT-30MLA: 243 BCT-30MA: 222 BCT-20MA: 150 BCT-10MA: 78

BCT-5MA: 42

Weight (g)

BCT-90MLA: 790
 BCT-60MLA: 730
 BCT-30MLA: 660
 BCT-30MA: 276
 BCT-20MA: 256
 BCT-10MA: 241
 BCT-5MA: 231

Magnetic Properties

Intrinsic coercivity (kA/m)	115.4
-----------------------------	-------

Retentivity (mT)	225.0
------------------	-------

Squareness (Br/Bm)	0.78
--------------------	------

Physical Properties

Thickness: backcoating (μm)	1.0
-----------------------------	-----

Thickness: base (μm)	10.0
----------------------	------

Thickness: magnetic layer (μm)	3.5
--------------------------------	-----

Thickness: total (μm)	14.5
-----------------------	------

Width (mm)	12.65
------------	-------

Gallery

