

CU-F31

Dynamic Super-cardioid Mic Capsule for use with DWM-02 and DWZ series handheld transmitters



Overview

Designed for vocal applications including live music performance, the CU-F31 is a dynamic super-cardioid mic capsule designed for use with the DWM-02 and DWZ series handheld transmitters

Dynamic Microphone Capsule for use with DWM-02 Handheld Transmitter

Specifications

General Specifications	
Capsule Type	Dynamic
Dimensions	1 7/8" Diameter x 3 1/8" Height
Directivity	Super-cardioid
Frequency Response	60 Hz - 18 kHz
Weight	5.1 oz.

Gallery

