

The smart digital imaging platform for medical environments



Overview

Manage, share and store HD, 4K, and 3D* video content across hospital networks and improve workflows

Features

Seamless, intuitive image workflow

Scalable, modality neutral and future-proof

Compatible with medical image sources

User-focused video editing

Flexible display options

Scalable NUCLeUS Smart Applications

Teaching, training and knowledge sharing

Streamlined information management

Flexible integration and remote support

Open for developers

Reducing patient anxiety in the OR

Protecting patient privacy

Related products



LMD-X310MT

31-inch 4K 3D/2D LCD medical monitor



LMD-X550MT

55-inch 4K 3D/2D LCD medical monitor



LMD-X310MD

31-inch 4K 2D LCD medical monitor



LMD-X550MD

55-inch 4K 2D LCD medical monitor



FW-55BZ35F

55" BRAVIA 4K HDR Professional Display



SRG-X120

IP 4K* Pan-Tilt-Zoom Camera with NDI®**|HX capability



SRG-X400

IP 4K* Pan-Tilt-Zoom Camera with NDI®**|HX capability



FW-75BZ35F

75" BRAVIA 4K HDR Professional Display



LMD-X2705MD

27-inch 4K 2D LCD medical monitor

Gallery

