

## UPC-24LA

Self-laminating colour printing pack (L-size) for UP-25MD and UP-D25MD medical printers



### Overview

A6 (L-size) self-laminating dye sublimation colour printing pack for use with Sony UP-25MD and UP-D25MD digital colour printers.

Pack contains 160 sheets of print media and ribbon (40 sheets x 4 packs).

Media size 144 x 100 mm.

### Specifications

#### Specifications

Type of Media	Self-Laminating Color Printing Pack
Contents	160 sheets of print paper 4 roll of ink ribbon
Paper/Film Size	144 x 100 mm (A6) 5 5/8 x 4 inches
Print Pages for One	160 sheets of print paper (40sheets x 4)

Roll

4 roll of ink ribbon

---

Supported Printers

UP-25MD, UP-D25MD

---

## Gallery

