

## MB-L22

Mounting bracket for LMD-B240/LMD-A220/PVM-A250 monitors



### Overview

**Mounting bracket for LMD-B240 24-inch basic grade LCD monitor, LMD-A220 22-inch high grade LCD monitor or PVM-A250 25-inch TRIMASTER EL™ OLED high grade picture monitor**

The MB-L22 is a mounting bracket specifically designed for mounting the LMD-B240/LMD-A220/PVM-A250 Sony monitor directly on a 19-inch EIA standard rack.

#### **Secure mounting in 19-inch EIA standard rack**

Mounting bracket for easy monitoring in operations room and OB truck.

### Features

#### **Metal mounting bracket**

Designed to attach to LMD-B240, LMD-A220 or PVM-A250 monitor to mount directly in a 19-inch EIA standard rack.

PVM-A250: Height 10U

LMD-B240, LMD-A220: Height 9U

## Gallery

