

## UPC-510

Colour printing pack for UP-50, UP-D50 and UP-51MD medical printers



### Overview

A5 dye sublimation colour printing pack for use with Sony UP-50, UP-D50 and UP-51MD printers.

Pack contains 200 sheets of print media plus ink ribbon.

Media size 148 x 210 mm.

### Specifications

#### Specifications

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| Type of Media            | Color Printing Pack                               |
| Contents                 | 200 sheets of print paper<br>1 roll of ink ribbon |
| Paper/Film Size          | 210mm x 148 mm (A5)<br>8 3/8 x 5 7/8 inches       |
| Print Pages for One Roll | 200 sheets of print paper<br>1 roll of ink ribbon |
| Supported Printers       | UP-51MD,UP-51MDP,UP-D50                           |



## Gallery

