

CCMC-T05

5 metre camera cable for MCC-3000MT



Overview

Extends flexibility of medical camera

The CCMC-T05 camera cable allows you to extend the MCC-3000MT medical two-piece camera. The MCC-3000MT has a compact camera head unit, that can be attached easily either to a surgical microscope, a surgical boom arm, or to a ceiling, and a separated camera control unit. The 5 m cable increases the distance between the camera head and the camera control unit.

MDD approved

The camera cable meets MDD approval in the EU.

Features

Extends distance between camera head and camera control unit

The 5 m camera cable is fully MMD approved in the EU and extends the flexibility of the MCC-3000MT in surgical conditions.

Gallery

