

## MDR-7510

Studio professional headphones



Everything Audio  
Network  
STELLAR SOUND!

### Overview

#### **Precision made for precision work**

An excellent choice for both reference monitoring and audio mixing applications, the 7510 extends the MDR-7500 line with a new fixed design of rugged construction. With 50mm drivers utilizing PET (Polyethylene Terephthalate) diaphragms, these headphones provide a wide frequency response, broad dynamic range, and accurate sound reproduction.

The 7510 offers high power handling capability and uses neodymium magnets and OFC (oxygen-free copper) voice coils. Noise isolation earpads with pliable cushion material conform to the natural shape of the ear and provide comfort and acoustic isolation.

### Features

#### **50mm Driver Unit with PET Diaphragm**

#### **Closed-Ear Design**

#### **Wide Frequency Response**

#### **High Power Handling Capacity 2,000mW**

#### **Stereo Unimatch Plug**

#### **Gold Connectors and OFC Cord**

**Supplied Soft Case**

## Specifications

| Group heading      |  |
|--------------------|--|
| Headphone Type     | Dynamic, Closed                              |
| Magnet Type        | Neodymium                                    |
| Driver Size        | 50.0mm                                       |
| Frequency Response | 5-40kHz                                      |
| Impedance          | 24 Ohms                                      |
| Sensitivity        | 108 dB/mW                                    |
| Power Handling     | 2,000mW                                      |
| Plug Type          | Gold, Stereo Unimatch<br>3.5mm/6.3mm adapter |
| Cord Length        | 3m (9.8 ft)                                  |
| Weight             | 261g (9.2 oz)                                |

## Gallery

