

## MDR-7550

High quality in-ear monitoring headphones for live stage use



### Overview

#### **Precision made for precision work**

Ideally suited for musicians who listen to custom-crafted mixes of vocals and stage instrumentation during live performance, this model utilizes large calibre 16mm dynamic drivers to enable a wide dynamic range with clear and accurate sound reproduction.

The 7550 offers a unique multi-layer diaphragm (comprised of hundreds of layers of polymer, each less than 0.01-micron thick) and a flexible ear hanger made of TEKNOROTE™ material, which easily folds and maintains its shape, insuring a superior fit for each individual ear. Also included are newly developed hybrid silicone rubber earbuds, with three sizes to fit different ear types.

### Features

#### **16mm Driver Unit with ML Diaphragm**

#### **In-Ear Monitor (IEM) Design**

#### **Flexible Ear Hanger Insures Superior Fit**

#### **Hybrid Silicone/Rubber Earbuds, 3 Sizes Supplied**

#### **Supplied Soft Case**

## Specifications

## Group heading

Headphone Type	Dynamic, In-Ear
Magnet Type	Neodymium
Driver Size	16.0 mm
Frequency Response	3-28kHz
Impedance	16 Ohms
Sensitivity	108 dB/mW
Power Handling	500mW
Plug Type	Gold, Stereo Unimatch 3.5mm/6.3mm adapter
Cord Length	1.6m (5.25 ft)
Weight	7.09g (.25 oz)

## Gallery

