

HVO-AUDIOKIT

Audio narration kit for the
HVO500/550MD Medical Grade
DVD Recorder



Overview

Audio narration kit for the HVO500/550MD Medical Grade DVD Recorder

This kit will give the HVO500/550MD an audio feature enabling an audio track to be recorded onto a DVD as video is being recorded. This feature can later be used for review and training purposes.

Highlights

- Rane Microphone Amplifier
- Sony ECM-673 Microphone
- XLR/RCS Adapter
- XLR male to XLR female cable
- RCA/RCA Y Cable

Gallery

