

## XDBK-102

Optional 250 GB SSD Drive for RAID-4 Configurations (XDS-PD2000 only)



### Overview

**Third 250 GB solid state drive in a drive carrier to configure XDS-PD2000 internal storage as a RAID-4 configuration**

The XDBK-102 is a 250 GB solid state drive with drive sled that enables you to further enhance the reliability of the Professional Media Station's internal storage by creating a RAID-4 configuration.

**Reliable storage solution for XDS-PD2000**

### Features

**Internal storage for the XDS-PD2000 Media Station**

The FSR 250 GB SSD Drive for RAID-4 Configuration further enhances the reliability of your XDS-PD2000 Professional Media Station's internal storage in a RAID-4 set-up.

## Gallery

