

LMD-SN10

Network Management software
for LMD-DM series Medical
Diagnostic Displays



Overview

This easy-to-use PC software allows monitoring of multiple LMD-DM series diagnostic displays via a hospital network and remote terminal.

Centralised management functions include: notice of local display malfunction; display profiles (model name, serial number, error code/log); workstation configuration in network; and maintenance/consistency test & report.

Compatible with LMD-DM50, LMD-DM30, LMD-DM20, LMD-DM30C, LMD-DM20C medical diagnostic displays.

Features

Wide range of management functions for networked diagnostic displays

Functions include: notice of local display malfunction; display profiles (model name, serial number, error code/log); workstation configuration in network; and maintenance/consistency test & report.

Includes 10 licenses

Allows networked management of up to ten LMD-DM series diagnostic displays.

Gallery

