

## UPA-500

Cleaning Kit for for UP-DF750,  
UP-DF550, UP-DF500 series  
medical printers



### Overview

This cleaning kit is compatible with Sony UP-DF750, UP-DF550 and UP-DF500 series medical printers.

Cleaning kit contents: 5 cleaning rollers, 5 sheets of cleaning paper, 1 sheet of wrapping film.

### Specifications

#### Specifications

Description	Cleaning Kit
Contents	Cleaning roller x 5 pcs Cleaning paper x 5 pcs Wrapping Film x 1 pcs
Supported Printers	UP-DF500, UP-DF550, UP-DF750

## Gallery

