

## UPA-500

Cleaning Kit for for UP-DF750,  
UP-DF550, UP-DF500 series  
medical printers



### Overview

This cleaning kit is compatible with Sony UP-DF750, UP-DF550 and UP-DF500 series medical printers.

Cleaning kit contents: 5 cleaning rollers, 5 sheets of cleaning paper, 1 sheet of wrapping film.

### Specifications

#### Specifications

|                    |  |
|--------------------|--|
| Description        | Cleaning Kit   |
| Contents           | Cleaning roller x 5 pcs<br>Cleaning paper x 5 pcs<br>Wrapping Film x 1 pcs |
| Supported Printers | UP-DF500, UP-DF550, UP-DF750   |

## Gallery

