

## CU-F31

Dynamic Super-cardioid Mic Capsule for use with DWM-02 and DWZ series handheld transmitters



### Overview

Designed for vocal applications including live music performance, the CU-F31 is a dynamic super-cardioid mic capsule designed for use with the DWM-02 and DWZ series handheld transmitters

### **Dynamic Microphone Capsule for use with DWM-02 Handheld Transmitter**

### Specifications

#### General Specifications

|                    |                                 |
|--------------------|---------------------------------|
| Capsule Type       | Dynamic                         |
| Dimensions         | 1 7/8" Diameter x 3 1/8" Height |
| Directivity        | Super-cardioid                  |
| Frequency Response | 60 Hz - 18 kHz                  |
| Weight             | 5.1 oz.                         |

## Gallery

