

UPC-21S

A6 Small-size colour print pack for use with UP-27MD, UP-D25MD, UP-25MD, UP-D23MD, UP-21MD, UP-D21MD, UP-20 printers.



Overview

Small-size "A6" colour print pack.

This product is available in multiples of 6 units (1 master carton).

Mainly used in Endoscopy and Ultrasound applications.

Specifications

Generic Specifications

Pack contents	240 sheets of print media, 3 colour ink ribbons
---------------	---

Order unit	1 pack (multiples of 6 units)
------------	-------------------------------

Media size	100 x 90 mm
------------	-------------

Gallery

