

RM-91

Wired Remote Commander for medical printers



Overview

This optional wired remote commander can be used to control Sony UP-991AD, UP-971AD, UP-990AD, UP-970AD, UP-X898MD, UP-897MD and UP-25MD medical printers.

Connector: Stereo mini jack

Cable length: 5m


Weight: 80g

Specifications

Specifications

Description	Remote Commander
Connector	Stereo mini
Cable Length	5 m 16 feet 5 inches
Mass	80 g 2 7/8 oz
Supplied Accessories	Operation manual (1)
Supported Printers	UP-20, UP-21MD, UP-25MD, UP-55MD, UP-897MD, UP-990AD, UP-970AD

Gallery



dimension are approximate. All trademarks are the property of their respective owners.