

LMD-KT10

Calibration Kit for LMD-DM series
Medical Diagnostic Displays



Overview

This easy-to-use PC software and sensor kit supports consistent display quality control with Sony LMD-DM series monitors – an essential requirement with medical diagnostic display applications.

Compatible with LMD-DM50, LMD-DM30, LMD-DM20, LMD-DM30C, LMD-DM20C medical diagnostic displays.

Features

Display Calibrator

Allows easy, accurate display calibration with supplied calibration sensor.

Allows easy, accurate display calibration with supplied calibration sensor.

Allows accurate measurement of display luminance.

Display Quality Controller

Provides acceptance and consistency tests according to AAPM, DIN and JESRA quality assurance guidelines.

Display Utility Software

Allows central control and management of display information by remote system administrator.

Gallery

