

## LMD-SN10

Network Management software  
for LMD-DM series Medical  
Diagnostic Displays



### Overview

This easy-to-use PC software allows monitoring of multiple LMD-DM series diagnostic displays via a hospital network and remote terminal.

Centralised management functions include: notice of local display malfunction; display profiles (model name, serial number, error code/log); workstation configuration in network; and maintenance/consistency test & report.

Compatible with LMD-DM50, LMD-DM30, LMD-DM20, LMD-DM30C, LMD-DM20C medical diagnostic displays.

### Features

#### **Wide range of management functions for networked diagnostic displays**

Functions include: notice of local display malfunction; display profiles (model name, serial number, error code/log); workstation configuration in network; and maintenance/consistency test & report.

#### **Includes 10 licenses**

Allows networked management of up to ten LMD-DM series diagnostic displays.

## Gallery

