

## NSBK-HS05

2TB hard drive for the NSR-500



### Overview

#### **Provides 2T storage for the 16 channel Full HD network surveillance server**

Extends the capacity of the NSR-500 network surveillance recording server. With a 12 TB HDD (2 TB x 6), you can record more than 300 hours of high-quality full-HD images compressed at 4 Mbps with a frame rate of 30 fps.

#### **Long-duration recording**

Extends the capacity of the NSR-500 network surveillance recording server.

### Features

#### **Long-duration recording**

Extends the capacity of the NSR-500 network surveillance recording server. With a 12 TB HDD (2 TB x 6), you can record more than 300 hours of high-quality full-HD images compressed at 4 Mbps with a frame rate of 30 fps.

## Gallery

