

## UPC-21L

A6 Large-size colour print pack for use with UP-27MD, UP-D25MD, UP-25MD, UP-D23MD, UP-21MD, UP-20 printers.



### Overview

Large-size "A6" colour print pack.

This product is available in multiples of 6 units (1 master carton).

Mainly used in Endoscopy and Ultrasound applications.

### Specifications

#### Generic Specifications

Pack contents	200 sheets of print media, 4 colour ink ribbons
Order unit	1 pack (multiples of 6 units)
Media size	144 x 100 mm

## Gallery

