

# UPC-24SA

Self-laminating colour printing pack (S-size)  
for UP-25MD and UP-D25MD medical  
printers



## Overview

A6 (S-size) self-laminating dye sublimation colour printing pack for use with Sony UP-25MD and UP-D25MD digital colour printers.

Pack contains 180 sheets of print media and ribbon (60 sheets x 3 packs).

Media size 100 x 90 mm.

## Specifications

### Specifications

Type of Media	Self-Laminating Color Printing Pack
Contents	180 sheets of print paper 3 roll of ink ribbon
Paper/Film Size	100 x 90 mm (A7) 4 x 3 1/2 inches
Print Pages for One Roll	180 sheets of print paper (60sheets x 3) 3 roll of ink ribbon
Supported Printers	UP-25MD, UP-D25MD

Information on EU Data Act: [Click here](#) to check if your Sony product and/or related service is subject to the EU Data Act.

## Gallery



dimension are approximate. All trademarks are the property of their respective owners.