

## SU-560

Stand specifically designed for medical monitors, working with LMD-1951MD, LMD-2451MD and PVM-2551MD



### Overview

The Sony SU-560 is specifically designed for the LMD-1951MD, LMD-2451MD and PVM-2551MD medical grade LCD monitors. It features four threaded holes on the bottom to facilitate its permanent installation on a medical equipment cart.

### Features

#### **Stand specifically designed for medical monitors**

The SU-560 has been specifically designed for the LMD-1951MD, LMD-2451MD and PVM-2551MD medical grade LCD monitors.

#### **Convenient for equipment cart installation**

The SU-560 features four threaded holes on the bottom to facilitate its permanent installation on a medical equipment cart.

## Gallery

