

CU-C31

Condenser cardioid Mic Capsule
for use with DWM-02 and DWZ
series handheld transmitters



Overview

Designed for vocal applications including live music performance, the CU-C31 is a condenser cardioid mic capsule designed for use with the DWM-02 and DWZ series handheld transmitters

Condenser Microphone Capsule for use with DWM-02 Handheld Transmitter

Specifications

General Specifications

| | |
|--------------------|---------------------------------|
| Capsule Type | Condenser |
| Dimensions | 1 7/8" Diameter x 2 1/2" Height |
| Directivity | Cardioid |
| Frequency Response | 60 Hz - 20 kHz |
| Weight | 4.8 oz. |

Gallery

