

CU-F32

Dynamic Wide-Cardioid Mic Capsule for use with DWM-02 and DWZ series handheld transmitters



Overview

Designed for presentation applications, the CU-F32 is a dynamic wide-cardioid mic capsule designed for use with the DWM-02 handheld transmitter and DWZ series handheld transmitters

Dynamic Microphone Capsule for use with DWM-02 Handheld Transmitter

Specifications

General Specifications

| | |
|--------------------|---------------------------------|
| Capsule Type | Dynamic |
| Dimensions | 1 7/8" Diameter x 3 1/8" Height |
| Directivity | Wide-cardioid |
| Frequency Response | 70 Hz - 18 kHz |
| Weight | 5.1 oz. |

Gallery

