

## UPC-24SA

Self-laminating colour printing pack (S-size) for UP-25MD and UP-D25MD medical printers



### Overview

A6 (S-size) self-laminating dye sublimation colour printing pack for use with Sony UP-25MD and UP-D25MD digital colour printers.

Pack contains 180 sheets of print media and ribbon (60 sheets x 3 packs).

Media size 100 x 90 mm.

### Specifications

#### Specifications

|                     |   |
|---------------------|---|
| Type of Media       | Self-Laminating Color Printing Pack               |
| Contents            | 180 sheets of print paper<br>3 roll of ink ribbon |
| Paper/Film Size     | 100 x 90 mm (A7)<br>4 x 3 1/2 inches              |
| Print Pages for One | 180 sheets of print paper (60sheets x 3)          |

Roll

3 roll of ink ribbon

---

Supported Printers

UP-25MD, UP-D25MD

---

## Gallery

