

## CU-F32

Dynamic Wide-Cardioid Mic Capsule for use with DWM-02 and DWZ series handheld transmitters



### Overview

Designed for presentation applications, the CU-F32 is a dynamic wide-cardioid mic capsule designed for use with the DWM-02 handheld transmitter and DWZ series handheld transmitters

### **Dynamic Microphone Capsule for use with DWM-02 Handheld Transmitter**

### Specifications

#### General Specifications

Capsule Type	Dynamic
Dimensions	1 7/8" Diameter x 3 1/8" Height
Directivity	Wide-cardioid
Frequency Response	70 Hz - 18 kHz
Weight	5.1 oz.

## Gallery

