

## UPC-21L

A6 Large-size colour print pack for use with UP-27MD, UP-D25MD, UP-25MD, UP-D23MD, UP-21MD, UP-20 printers.



### Overview

Large-size "A6" colour print pack.

This product is available in multiples of 6 units (1 master carton).

Mainly used in Endoscopy and Ultrasound applications.

### Specifications

#### Generic Specifications

|               |   |
|---------------|---|
| Pack contents | 200 sheets of print media, 4 colour ink ribbons |
| Order unit    | 1 pack (multiples of 6 units)                   |
| Media size    | 144 x 100 mm                                    |

## Gallery

