

UPC-24LA

Self-laminating colour printing pack (L-size) for UP-25MD and UP-D25MD medical printers



Overview

A6 (L-size) self-laminating dye sublimation colour printing pack for use with Sony UP-25MD and UP-D25MD digital colour printers.

Pack contains 160 sheets of print media and ribbon (40 sheets x 4 packs).

Media size 144 x 100 mm.

Specifications

Specifications

| | |
|---------------------|---|
| Type of Media | Self-Laminating Color Printing Pack |
| Contents | 160 sheets of print paper 4 roll of ink ribbon |
| Paper/Film Size | 144 x 100 mm (A6) 5 5/8 x 4 inches |
| Print Pages for One | 160 sheets of print paper (40sheets x 4) |

Roll

4 roll of ink ribbon

Supported Printers

UP-25MD, UP-D25MD

Gallery

