

## UPC-510

Colour printing pack for UP-50, UP-D50 and UP-51MD medical printers



### Overview

A5 dye sublimation colour printing pack for use with Sony UP-50, UP-D50 and UP-51MD printers.

Pack contains 200 sheets of print media plus ink ribbon.

Media size 148 x 210 mm.

### Specifications

#### Specifications

Type of Media	Color Printing Pack
Contents	200 sheets of print paper 1 roll of ink ribbon
Paper/Film Size	210mm x 148 mm (A5) 8 3/8 x 5 7/8 inches
Print Pages for One Roll	200 sheets of print paper 1 roll of ink ribbon
Supported Printers	UP-51MD,UP-51MDP,UP-D50



## Gallery

