

MDR-7510

Studio professional
headphones



Overview

An excellent choice for both reference monitoring and audio mixing applications, the 7510 extends the MDR-7500 line with a new fixed design of rugged construction. With 50mm drivers utilizing PET (Polyethylene Terephthalate) diaphragms, these headphones provide a wide frequency response, broad dynamic range, and accurate sound reproduction.

The 7510 offers high power handling capability and uses neodymium magnets and OFC (oxygen-free copper) voice coils. Noise isolation earpads with pliable cushion material conform to the natural shape of the ear and provide comfort and acoustic isolation.

Features

Specifications

Group heading

Headphone Type	Dynamic, Closed
Magnet Type	Neodymium
Driver Size	50.0mm
Frequency Response	5-40kHz
Impedance	24 Ohms
Sensitivity	108 dB/mW
Power Handling	2,000mW
Plug Type	Gold, Stereo Unimatch 3.5mm/6.3mm adapter
Cord Length	3m (9.8 ft)
Weight	261g (9.2 oz)

Gallery

