

CFV-50SC

Shield Cover



Overview

This pack contains 50 x disposable hygienic single-use covers (cover only, eye shield and shield frame are not included). Each cover is intended to be worn in conjunction with the Sony CFV-E30SK 3D eye shield.

Enabling unimpaired 3D viewing in the Operating Room, the transparent disposable cover protects the eye shield from dirt, blood and other fluid splashes while viewing images on compatible Sony 3D medical displays.

The CFV-50SC is compatible with the CFV-E30SK(3D), CFV-E20SK(2D) eye shields and these Sony 3D medical displays: LMD-X310MT, LMD-X550MT.

Please note: This product is not intended for use as personal protective equipment.

Specifications

General

| | |
|-----------------------|---------------|
| Operating Temperature | 0 °C to 40 °C |
|-----------------------|---------------|

| | |
|--------------------|------------|
| Operating Humidity | 30% to 85% |
|--------------------|------------|

| | |
|--------------------------------------|---|
| Storage/Transport Temperature | -20°C to +60°C |
| Storage/Transport Humidity | 20% to 90% |
| Operating/Storage/Transport Pressure | 700 hPa to 1,060 hPa |
| Dimensions *1 | 250 x 204 mm |
| Mass | 8.5 g (1) 425 g (50) 0.299 oz (1) 14.991 oz (50) |
| Supplied Accessories | Operating Instructions (1) |

*1 The values for dimensions are approximate.

Related products



CFV-E30SK
3D Eye Shield Kit



LMD-X310MT
31-inch 4K 3D/2D LCD medical monitor



LMD-X550MT
55-inch 4K 3D/2D LCD medical monitor

Gallery

