

SU-560

Stand specifically designed for medical monitors, working with LMD-1951MD, LMD-2451MD and PVM-2551MD



Overview

The Sony SU-560 is specifically designed for the LMD-1951MD, LMD-2451MD and PVM-2551MD medical grade LCD monitors. It features four threaded holes on the bottom to facilitate its permanent installation on a medical equipment cart.

Features

Stand specifically designed for medical monitors

The SU-560 has been specifically designed for the LMD-1951MD, LMD-2451MD and PVM-2551MD medical grade LCD monitors.

Convenient for equipment cart installation

The SU-560 features four threaded holes on the bottom to facilitate its permanent installation on a medical equipment cart.

Gallery

