

AC-120MD

AC Adapter for Medical Monitors



Overview

Provides DC power for compatible Sony medical monitors

The AC-120MD is an AC adapter that's designed for use with these compatible Sony LCD medical monitors only:

LMD-2735MD/2735MC (AC-120MD available as optional accessory)

LMD-2435MD/2435MC (AC-120MD available as optional accessory)

LMD-2765MD/2765MC (AC-120MD included as standard with monitor)

LMD-2760MD/2760MC (AC-120MD included as standard with monitor)

The adapter can operate wherever an AC mains outlet (100-240V, 50/60Hz) is available. It provides a regulated 24V DC output suitable for powering the monitors as detailed above.

Specifications

Specifications

Power

AC IN: 100 V-240 V, 50/60 Hz, 2.5 A-1.0 A

| | |
|--------------------------------------|--|
| | DC OUT: 24 V 6.25 A |
| Operating Temperature | 0 °C to 35 °C (32 °F to 95 °F) |
| Operating Humidity | 30% to 85% (no condensation) |
| Storage/Transport Temperature | -20 °C to +60 °C (-4 °F to +140 °F) |
| Storage/Transport Humidity | 20% to 90% |
| Operating/Storage/Transport Pressure | 700 hPa to 1060 hPa |
| Dimensions (W x H x D) *1 | 170 × 40 × 85 mm 6 3/4 × 1 5/8 × 3 3/8 inches |
| Mass | 1.2 kg 2 lb 10 oz |
| Supplied Accessories | Instructions for Use (1) Service Contact List (1) |

Note

*1

The values for dimensions are approximate.

Gallery

