

## CCMC-T05

5 metre camera cable for MCC-3000MT



### Overview

#### **Extends flexibility of medical camera**

The CCMC-T05 camera cable allows you to extend the MCC-3000MT medical two-piece camera. The MCC-3000MT has a compact camera head unit, that can be attached easily either to a surgical microscope, a surgical boom arm, or to a ceiling, and a separated camera control unit. The 5 m cable increases the distance between the camera head and the camera control unit.

#### **MDD approved**

The camera cable meets MDD approval in the EU.

### Features

#### **Extends distance between camera head and camera control unit**

The 5 m camera cable is fully MMD approved in the EU and extends the flexibility of the MCC-3000MT in surgical conditions.

## Gallery

