

## UPC-R81MD

Self-laminating colour printing pack (US letter size) for UP-DR80MD medical printer



### Overview

US letter size dye sublimation colour printing pack for use with Sony UP-DR80MD digital colour printer.

Self-laminating, curl free media for high print durability and longevity. Protective laminated coating is resistant to light, heat, water, humidity and fingerprints.

1 roll produces 50 prints.

1 pack contains 2 ribbons and 2 rolls of media, to yield 100 prints per pack.

Media length per roll 216 mm x 15 m (8 1/2 x 590 inches).

### Specifications

#### Specifications

Type of Media	Self-Laminating Color Printing Pack
---------------	-------------------------------------

Contents	Paper roll for 50 prints x 2
----------	------------------------------

Paper/Film Size	216 mm x 15 m 8 1/2 x 590 inches
-----------------	-------------------------------------

---

Print Pages for One  
Roll

50 pages x 2/pack

---

Supported Printers

UP-DR80MD

---

## Related products



### **UP- DR80MD**

A4 Digital Colour  
Printer

## Gallery

