

## CU-C31

Condenser cardioid Mic Capsule  
for use with DWM-02 and DWZ  
series handheld transmitters



### Overview

Designed for vocal applications including live music performance, the CU-C31 is a condenser cardioid mic capsule designed for use with the DWM-02 and DWZ series handheld transmitters

### **Condenser Microphone Capsule for use with DWM-02**

### **Handheld Transmitter**

### Specifications

#### General Specifications

Capsule Type	Condenser
Dimensions	1 7/8" Diameter x 2 1/2" Height
Directivity	Cardioid
Frequency Response	60 Hz - 20 kHz
Weight	4.8 oz.

## Gallery

