

## CU-F32

Dynamic Wide-Cardioid Mic Capsule for use with DWM-02 and DWZ series handheld transmitters



### Overview

Designed for presentation applications, the CU-F32 is a dynamic wide-cardioid mic capsule designed for use with the DWM-02 handheld transmitter and DWZ series handheld transmitters

### **Dynamic Microphone Capsule for use with DWM-02 Handheld Transmitter**

### Specifications

#### General Specifications

|                    |                                 |
|--------------------|---------------------------------|
| Capsule Type       | Dynamic                         |
| Dimensions         | 1 7/8" Diameter x 3 1/8" Height |
| Directivity        | Wide-cardioid                   |
| Frequency Response | 70 Hz - 18 kHz                  |
| Weight             | 5.1 oz.                         |

## Gallery

