

## UPT-210BL

Blue thermal transparent film for UP-990AD medical printer



### Overview

Blue thermal transparent film for use with Sony UP-991AD, UP-990AD and UP-980AD hybrid graphic printers.

Approximately 42 prints per roll.

Media size 210 mm (W) x 12.5 m.

### Specifications

#### Specifications

Type of Media	Blue Thermal Transparent Film
Contents	210 mm x 12.5 m
Paper/Film Size	210 mm x 12.5m 8 1/4 x 492 inches
Print Pages for One Roll	Approximately 42 prints
Supported Printers	UP-991AD, UP-990AD, UP-980AD

## Gallery

