

## CU-F32



Designed for presentation applications, the CU-F32 is a dynamic wide-cardioid mic capsule designed for use with the DWM-02 handheld transmitter and DWZ series handheld transmitters

**Dynamic Microphone Capsule for use with DWM-02  
Handheld Transmitter**

### General Specifications

Capsule Type	Dynamic
Dimensions	1 7/8" Diameter x 3 1/8" Height
Directivity	Wide-cardioid
Frequency Response	70 Hz - 18 kHz
Weight	5.1 oz.

