

UPP-210SE

Black and white thermal print media for UP-991AD / 990AD / 980AD / 971AD / 970AD / 960AD / 930AD / 910AD medical printers



Overview

Standard monochrome print media for use with Sony UP-991AD / 990AD / 980AD / 971AD / 970AD / 960AD / 930AD / 910AD Black and White Printers.

Approximately 130 prints per roll. (Picture size:150x210mm)

Media size 210 mm (W) x 25 m.

Specifications

| Specifications | |
|--------------------------|--|
| Type of Media | High Quality Printing Paper |
| Contents | 210 mm x 25 m |
| Paper/Film Size | 210 mm x 25 m 8 1/4 x 984 inches |
| Print Pages for One Roll | Approximately 130 prints per roll (Picture size:150x210mm) |
| Supported Printers | UP-990AD, UP-971AD, UP-970AD |

Gallery

